


OUTSTANDING HYGIENE

Also suitable
for
AWT systems.




Sterile equipment transport trolleys, open



| | Description | Dimensions (WxDxH) | Art. No. |
|--|---|----------------------|----------|
| Open sterile equipment transport trolley with shelves* | | | |
| | Open 25 x 25 mm square tube frame construction made of stainless steel With 3 height-adjustable shelves (max. loading capacity 60 kg per shelf) 4 electrically conductive steering castors, 2 of which have brakes, 125 mm dia., with wall guards | | |
| 1 | For holding a maximum of 6 STU sterile products (for a maximum of 3 small shelves of 679 x 564 x 30 mm each) | 760 x 695 x 1397 mm | 543 137 |
| 2 | For holding a maximum of 9 STU sterile products (for a maximum of 3 large shelves of 1009 x 564 x 30 mm each) | 1090 x 695 x 1397 mm | 543 139 |
| Open sterile equipment transport trolley with support brackets* | | | |
| | Open 25 x 25 mm square tube frame construction made of stainless steel With 6 height-adjustable support brackets (max. loading capacity 20 kg per pair of support brackets) 4 electrically conductive steering castors, 2 of which have brakes, 125 mm dia., with wall guards | | |
| 3 | For holding a maximum of 6 STU sterile products with 6 pairs of support brackets (incl. 1 STW frame for support brackets) | 760 x 695 x 1397 mm | 543 138 |
| 4 | For holding a maximum of 9 STU sterile products with 9 pairs of support brackets (incl. 2 STW frames for support brackets) | 1090 x 695 x 1397 mm | 543 140 |

* Can be easily modified in just seconds – to shelves or support brackets depending on the model – completely without tools

Sterile equipment transport trolleys, closed



| | Description | Exterior dimensions (WxDxH) Inside dimensions (WxDxH) | Empty weight | Art. No. |
|---|---|--|--------------|----------------|
| Closed sterile equipment transport trolley with shelves* | | | | |
| Closed, single-walled model in stainless steel, 2 double-walled hinged door with all-round seals and stop bar, door in opening angle 270°, can be locked at 255°, with 2 removable, height-adjustable shelves, 2 vertical push handles on both sides, 4 corner guards, 2-point locking, 4 galvanised steering castors, 160 mm dia., 2 of which have brakes, loading weight 180 kg, max. loading capacity per shelf: 60 kg | | | | |
| 1 | For holding a maximum of 6 STU sterile products | 860 x 785 x 1485 mm 634 x 622 x 1132 mm | 94.5 kg | 543 131 |
| 2 | For holding a maximum of 9 STU sterile products | 1190 x 785 x 1485 mm 964 x 622 x 1132 mm | 117.5 kg | 543 133 |
| Closed sterile equipment transport trolley with support brackets* | | | | |
| Closed, single-walled model in stainless steel, 2 double-walled hinged door with all-round seals and stop bar, door in opening angle 270°, can be locked at 255°, 2 vertical push handles on both sides, 4 corner guards, 2-point locking, 4 galvanised steering castors, 160 mm dia., 2 of which have brakes, loading weight 180 kg, max. loading capacity per pair of support brackets: 20 kg | | | | |
| 3 | For holding a maximum of 6 STU sterile products with 6 pairs of height-adjustable support brackets (incl. 1 STW frame for support brackets, 2 connecting profiles for STW frame and 1 fixing profile for STW frame) | 860 x 785 x 1485 mm 634 x 622 x 1132 mm | 102.5 kg | 543 132 |
| 4 | For holding a maximum of 9 STU sterile products with 9 pairs of height-adjustable support brackets (incl. 2 STW frames for support brackets, 4 connecting profiles for STW frame and 2 fixing profiles for STW frame) | 1190 x 785 x 1485 mm 964 x 622 x 1132 mm | 130.5 kg | 543 134 |

Sterile equipment transport trolleys, closed and fully welded

| | Description | Exterior dimensions (WxDxH) Inside dimensions (WxDxH) | Empty weight | Art. No. |
|--|--|--|--------------|----------------|
| Closed sterile equipment transport trolley, fully welded, with shelves* | | | | |
| Closed, fully-welded model in stainless steel, suitable for washing systems, 2 double-walled hinged door with all-round seals and stop bar, door in opening angle 270°, can be locked at 255°, with 2 removable, height-adjustable shelves, 2 vertical push handles on both sides, 4 corner guards, 2-point locking, 4 galvanised, electrically conductive steering castors, 160 mm, dia., 2 of which have brakes, for indoor/outdoor transport, loading weight 180 kg, max. loading capacity per shelf: 60 kg | | | | |
| | For holding a maximum of 6 STU sterile products | 870 x 780 x 1495 mm 634 x 622 x 1132 mm | 107 kg | 543 194 |
| | For holding a maximum of 9 STU sterile products | 1200 x 780 x 1495 mm 964 x 622 x 1132 mm | 113 kg | 543 196 |
| Closed sterile equipment transport trolley, fully welded, with support brackets* | | | | |
| Closed, fully-welded model in stainless steel, suitable for washing systems, 2 double-walled hinged doors with all-round seals and stop bar, door in opening angle 270°, can be locked at 255°, 2-point locking, 4 galvanised electrically conductive steering castors, 160 mm dia., 2 of which have brakes, for indoor/outdoor transport, loading weight 180 kg, max. loading capacity per pair of support brackets: 20 kg | | | | |
| | For holding a maximum of 6 STU sterile products with 6 pairs of height-adjustable support brackets (incl. 1 STW frame for support brackets, 2 connecting profiles for STW frame) | 870 x 780 x 1495 mm 634 x 622 x 1132 mm | 132 kg | 543 195 |
| | For holding a maximum of 9 STU sterile products with 9 pairs of height-adjustable support brackets (incl. 2 STW frame for support brackets, 4 connecting profiles for STW frame) | 1200 x 780 x 1495 mm 964 x 622 x 1132 mm | 144 kg | 543 197 |

* Can be easily modified in just seconds – to shelves or support brackets depending on the model – completely without tools